

Atlas Toilevator®

Features

- Adds 3.5" to toilet seat height
- Rated for users up to 500lb
- Discreet: Blends in with the toilet base
- Sanitary!
- Easy to install
- Kit contains all parts needed for installation
- Color: Paintable white



Technical specification

Atlas Toilevator®

Toilet Riser – Fits most round and elongated toilets

- 500lb (227kg) maximum weight capacity
- Adds 3.5" to the height of the base of the toilet
- Discreet – no one will realize it is there
- Fits toilets with the base as close as 7" (18cm) from the wall
- May also be painted to match your existing toilet, we suggest using a primer suitable for ABS plastic
- Will fit toilets where the center of the soil pipe is as close as 7" to the wall

| | |
|---------------------------------|--------------------|
| Product Code | ATV300 |
| Dimensions | 23.25"x11.44"x3.5" |
| Addition to the toilet's height | 3.5" |
| Weight Capacity* | white |
| Available Colors | 500lb |



Atlas Toilevator Grande®

Toilet Riser – Fits newer, larger base toilets

- 500lb (227kg) maximum weight capacity
- Adds 3.5" to the height of the base of the toilet
- Discreet – no one will realize it is there
- Fits toilets where the center of the soil pipe is as close as 7" to the wall
- Toilevator® can also be painted to match your existing toilet, we suggest using a primer suitable

| | |
|---------------------------------|--------------------|
| Product Code | ATV350 |
| Dimensions | 28.75"x15.25"x3.5" |
| Addition to the toilet's height | 3.5" |
| Weight Capacity* | white |
| Available Colors | 500lb |

