

Hercules Fixed Support and Therapy System



Features

- End-to-end SafetyGrip® on all vertical poles
- Horizontal tracks screw into wall studs for maximum strength & stability
- Rated for all users including bariatric up to 500lb
- Multiple options to make life easy
- Customize the installation that's perfect for YOU!

Technical specification

Hercules Circa Wall Mounted Stability Poles

- Ideal for around a bed or in a physical therapy room
- When paired with a Hygeia kit, facilitates exercise to gain strength, improve balance, and prevent muscle atrophy
- Can be installed on wooden or metal stud walls.
- Requires a minimum of 68" of unobstructed wall space, 60" to 82" from the floor

Hercules Circa KIT 1

The Hercules fixed support and therapy system

Kit Contents:

1 x RCSG501 – Circa Track	Product code	RCSG500
1 x RCSG502 – Circa Pole	Weight Capacity*	500lb
1 x RCSG503 – Circa Arm	Available Colors	white

Hercules Circa KIT 2

The Hercules fixed support and therapy system

Kit Contents:

1 x RCSG501 – Circa Track	Product code	RCSG510
2 x RCSG502 – Circa Pole	Weight Capacity*	500lb
2 x RCSG503 – Circa Arm	Available Colors	white

Hercules Circa Parts

	Hercules Circa Track	Hercules Circa Pole	Hercules Circa Arm
Product code	RCSG501	RCSG502	RCSG503
Pole Height	4"	4"	4"
Contents	2 x Track pieces	1 x Pole piece	2 x Arm pieces
Weight Capacity*	500lb	500lb	500lb
Available Colors	white	white	white

Hercules Hygeia



Hygeia Kit

Product code	HHSGA575
Available Colors	Black grip and straps 4 colored resistance bands included in each kit
Compatibility	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650 All Titan poles
	The Hygeia Exercise Kit Includes: 8 x Resistance bands (4 different strengths) and 3 sets of hand/ankle straps



Hygeia Clamp Kit

Product code	THSGA605
Available Colors	black, white
Compatibility	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650 All Titan poles
	Exercise Clamps Kit includes: 3 x Clamps for fastening therapy equipment

