

Hypnos Bed Rails



Features

- End-to-end SafetyGrip® for the most secure hold or the soft-touch Grip of the Original M-Rail
- Rated for users up to 300lb
- Dignified unassisted entry and exit from your bed
- Affordable, configurable, attractive
- Easy installation

Technical specification

Hypnos Morpheus Bed rail featuring SafetyGrip®

Rectangle double bar system

- Featuring SafetyGrip® for a strong, no-slip hold
- Provides support when getting into or out of bed and when moving from a lying position to a sitting position
- Two height settings, with a crossbar to help prevent entrapment
- Straps to the box spring of all bed sizes from twin to king
- Easy, quick assembly
- All surfaces protected by AlphaSan®
- Made of powder coated steel for a long life and easy leaning weighing up to 350lb.

Product code	HMSG450
Weight Capacity*	300lb
Available Colors	black, white, glow



Hypnos M-Rail® Bed rail

Soft-grip bed rail for box springs beds

- The original, patented M-Rail® with Soft Touch grip
- Provides support when getting into or out of bed and when moving from a lying position to a sitting position
- Two height settings, with a crossbar to help prevent entrapment
- Made of powder-coated steel for a long life and easy cleaning
- Straps to the box spring of all bed sizes from twin to king
- Easy, quick assembly

Product code	HMR400
Weight Capacity*	300lb
Available Colors	black, white



Hypnos Apolo Bed rail featuring SafetyGrip®

Horizontal Bed Rail

- Package includes clamps for a variety of bed frame styles including angle iron, clinical style inverted angle iron, and hospital style round, square or rectangular frames
- Made of powder-coated steel for a long life and easy cleaning
- Featuring SafetyGrip® for a strong, no-slip hold
- For use in home, clinical, or hospital settings
- Rotates 360° and locks in place at every 90°
- Fits any size metal frame bed.
- All surfaces protected by AlphaSan®

Product code	HASG100
Weight Capacity*	300lb
Available Colors	white



Hypnos Arcorail Bed Rail

Two Loop Bed Rail

- Package includes clamps for a variety of bed frame styles including angle iron, clinical style inverted angle iron, and hospital style round, square or rectangular frames
- Made of powder-coated steel for a long life and easy cleaning
- Dozens of possible hand positions
- For use in home, clinical, or hospital settings
- Rotates 360° and locks in place at every 90°
- Fits any size metal bed frame

Product code	HAR250
Weight Capacity*	300lb
Available Colors	white

