

Kratos Grab Bars

Features

- SafetyGrip® for the most secure hold
- Rated for up to 350lb
- Affordable, configurable, attractive
- Customize an installation that's perfect for YOU!




Technical specification

Kratos Capto Straight Wall mounted Grab Bar



Grab Bar length	14"	24"	36"	46"
Product code	KCSG714	KCSG724	KCSG736	KCSG746
Weight Capacity*	350lb	350lb	350lb	350lb
Available Colors	black, white	black, white	black, white	black, white

Kratos Capto Hinged Wall mounted Grab Bar



Capto V

Grab Bar length	2x24"
Product code	KCSG750
Weight Capacity*	350lb
Available Colors	black, white




Capto 9

Grab Bar length	2x24"
Product code	KCSG750
Weight Capacity*	350lb
Available Colors	black, white

Kratos Bia Hinged rear wall grab bar


- Support arm that folds neatly out of the way when not in use.
- Typically mounted on one or both sides of the toilet.
- Innovative offset rail helps reduce wrist strain.
- Protected by antimicrobial AlphaSan®.



Grab Bar length	32"
Product code	KBSG775
Weight Capacity*	350lb
Available Colors	black, white

Kratos Poseidon Bath mounted hand rail

- Protected by antimicrobial AlphaSan®
- Fully adjustable to suit most bath tubs
- Not suitable for fiberglass surfaces.



Grab Bar length	12"
Product code	KPSG780
Weight Capacity*	350lb
Available Colors	black, white

Kratos Trapeze Bath mounted hand rail

- Ceiling mounted trapeze
- Protected by antimicrobial AlphaSan®
- Designed to assist with balance, sitting up, positioning and stability for persons weighing up to 350lb.



Product code	KHSG706
Weight Capacity*	350lb
Available Colors	black, white

