

Titan Stability Pole Range


Features

- Multiple uses, from balance assistance to strength training & stretching
- Pairs perfectly with Therabands® and other therapy bands
- Tension mounted for go-anywhere convenience
- Allows 360° of unrestricted movement
- The grip area benefits from a full 27 inches of SafetyGrip®
- Easy one-person installation, tooling included
- Standard and bariatric models rated for users up to 500lb





Technical specification

Titan Poles



	Titan Pole	TitanOne	TitanPC	TitanHD
Product code	TSG600	TSG610	TSG615	TSG650
Maximum Height	8'3"	8'3"	8'3"	8'3"
Area of SafetyGrip®	27"	27"	27"	27"
Weight Capacity*	300lb	350lb	350lb	500lb
Available Colors	black, white, gray	black, white, gray	black, white, gray	black, white, gray

Titan Links



	Titan Extender	Titan - Rotobar Rotating Handle
Product code	TSG604	TSG601
Available Colors	black, white, gray	black, white
Compatibility	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650	TitanOne TSG610 TitanPC TSG615 TitanHD TSG650
	Height Extender Takes 8'3" Titan Pole to between 9"-3" to a maximum 10"-3" 7'9" to 8'3" - No height extension tubes required. 8'3" to 8'9" - 1 x 6" short extension tube required. 8'9" to 9'3" - 1 x 12" long extension tube required. 9'3" to 9'9" - 1 x 6" short and 1 x 12" long extension tube required. 9'9" to 10'3" - All 3 extension tubes required.	Safety Grip® Pole Rotating Grab Bar

Titan Therapy



	Titan Hygeia Kit	Titan Hygeia Clamp Kit	Titan Trapeze
Product code	HHSGA575	THSGA605	TSGA602
Available Colors	Black grip and straps 4 colored resistance bands included in each kit	black, white	white
Compatibility	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650 All Hercules poles	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650 All Hercules poles	Titan TSG600 TitanOne TSG610 TitanPC TSG615 TitanHD TSG650
	The Titan Exercise Kit Includes: 8 x Resistance bands (4 different strengths) and 3 sets of hand/ankle straps	Exercise Clamps Kit includes: 3 x Clamps for fastening therapy equipment	Pole mounted, height-adjustable trapeze Note: maximum weight capacity of 350lb, even when used with the TitanHD

