

# Toilevator®

## Sometimes All You Need Is a Little Height

The Toilevator® discreetly raises your toilet's height by 3.5 inches, making it easier to use without compromising the visual appeal of your bathroom. It's your key to a more comfortable, sanitary bathroom experience, ensuring peace of mind as you age in place or recover from an injury.



### Compare the Difference Between a Toilevator® and an Ordinary Toilet Seat Riser

See Why Toilevator® Beats the Competition



Toilevator®



Standard Seat Riser

	Toilevator	Standard Seat Riser
Maintains hygiene and cleanliness while reducing the need for frequent cleaning by eliminating user contact	✓	✗
Does not gain odors and unsightly staining from regular toilet use	✓	✗
Fastens safely and securely in place, preventing any sliding or misalignment that could cause risky falls	✓	✗
Discreetly and seamlessly blends with your toilet base	✓	✗
Can be painted to match bathroom decor	✓	✗




### What Our Customers Are Saying

*This is a great product for the elderly. I installed it for my sister and it's a shame I didn't find this many years ago. Great quality and very easy to install.*

*So much nicer to elevate the existing toilet than to use a bedside commode over a toilet as they often do at care centers.*

*This was wonderful for me. I am tall and have difficulty lifting myself off the toilet and the added height allows me the extra boost to stand.*

*Solid, stable, fit perfectly and easy to install. That's about as good as you can get!*



## Technical specifications



### Toilevator®

- 500lb (227kg) maximum weight capacity
- Adds 3.5" to the height of the base of the toilet
- Discreet - no one will realize it is there
- May also be painted to match your existing toilet, we suggest using a primer suitable for ABS plastic
- Will fit toilets where the center of the soil pipe is as close as 7" to the wall

<b>Product code</b>	TV300	<b>Addition to the toilet's height</b>	3.5"
<b>Dimensions</b>	23.25" x 11.44" x 3.5"	<b>Available Colors</b>	white



### Toilevator Grande® Fits larger base toilets

- All the same great features as the original Toilevator®, but built wider and longer to accommodate elongated toilet bases.

<b>Product code</b>	TV350	<b>Addition to the toilet's height</b>	3.5"
<b>Dimensions</b>	28.75" x 15.25" x 3.5"	<b>Available Colors</b>	white

