



SafetyGrip®
By Hartmobility

Titan Stability Poles

Grab. Hold. Feel The Difference.

Walking around your home shouldn't be a risky undertaking. With Titan Stability Poles, you can rest assured you'll have all the support you need.



| Titan Stability Poles

Whether in an institutional setting or at home, this floor-to-ceiling assist pole improves independence, mobility, and safety, allowing users to live freely and confidently.

Product Details

- Equipped with patented SafetyGrip® technology, ensuring the strongest grip possible.
- Developed with KITE, a world leader in complex rehabilitation science.
- Easy to install. Some models require no drilling and no fasteners.
- May be placed anywhere that a grab bar is needed. No wall required.
- Adjusts from 7'9" to 8'3," and can be paired with height extenders to reach as high as 10'3."
- Customize your installation with accessories: therapy band clamps, Rotobar rotating handle, height extenders, and more.




| Titan Links

Links are an economical, easy to install, inexpensive solution when a handrail is needed

Product Details

- Spans between two Titan poles or between a Titan pole and a wall bracket. No drilling or room modification necessary.
- May be used in sequence to achieve any desired horizontal handrail distance.
- Available in two lengths: 36" and 72."
- 36" hinged sections allow lift-up access for closets and doors.
- Works with all models of Titan poles.



| Titan Accessories

Customize your installation with accessories for added support, exercise and stretching assistance, or to accommodate the architecture of your home.

Product Details


- Roto-Bar attachment to extend support beyond standard stability pole and allow for easy repositioning.
- Trapeze system designed to assist with balance, transitioning to, and from, sitting, positioning, and stability.
- Resistance Band Clamps to aid exercise or rehabilitation. .



TITAN POLES

Color Code


- Black
- White
- Green



- Pole Diameter: 1.5"
- Grip Diameter: 1.75"
- Area of SafetyGrip®: 27"

Product Name	Titan Pole	TitanOne	TitanPC	TitanHD (Bariatric)	Height Extender Kit
					Takes 8'3" Titan Pole to between 9'3" to a maximum 10'3"
Weight Capacity*	300lb	350lb	350lb	500lb	
	<ul style="list-style-type: none"> • B • W • G 	<ul style="list-style-type: none"> • B • W • G 	<ul style="list-style-type: none"> • B • W • G 	<ul style="list-style-type: none"> • W 	<ul style="list-style-type: none"> • W • B
Product Code	TSG600	TSG610	TSG615	TSG650	TSG604

ACCESSORIES

	Product Name	Compatibility	Color	Product Code
	<p>Rotobar</p> <p>The Rotobar is an oval SafetyGrip® coated handle that utilizes pivot and lock technology to allow easy repositioning in 45 degree increments around the Titan Pole. The height placement of the Rotobar is adjustable, allowing it to fit most standard height requirements for the end-user</p>	<ul style="list-style-type: none"> • TSG610 • TSG615 • TSG650 	<ul style="list-style-type: none"> • W • B 	TSG601
	<p>Titan Hygeia Clamp Kit</p> <p>Turn your Titan pole into a home exercise and rehab tool with 3x clamps for fastening therapy bands</p>	<ul style="list-style-type: none"> • All Titan and Hercules poles 	<ul style="list-style-type: none"> • B • W 	THSGA605
	<p>Titan Trapeze</p> <p>Pole-mounted, height-adjustable trapeze Note: maximum weight capacity of 350lb, even when used with the TitanHD</p>	<ul style="list-style-type: none"> • All Titan and Hercules poles 	<ul style="list-style-type: none"> • W 	TSGA602

TITAN LINKS

	Product Name	Weight Capacity	Color	Product Code
	<p>Titan 72" Horizontal Rail</p> <p>Titan Link 72" long Horizontal Handrail - fixed.</p>	<ul style="list-style-type: none"> • 350lbs 	<ul style="list-style-type: none"> • W 	TLSG620
	<p>Titan 36" Horizontal Rail</p> <p>Titan Link 36" long Horizontal Handrail - fixed.</p>	<ul style="list-style-type: none"> • 350lbs 	<ul style="list-style-type: none"> • W 	TLSG625
	<p>Titan 36" Swing Gate</p> <p>The Titan Link 36" Swing Gate is a friction-assisted gate allowing access to doorways located on the opposite side of a Links handrail system.</p>	<ul style="list-style-type: none"> • 350lbs 	<ul style="list-style-type: none"> • W 	TSLG635
	<p>Titan Fixed Wall Mount Bracket</p> <p>The TLFM685 pairs with the 36" Swing Gate (TSLG635), allowing the Link System to terminate at a wall.</p>	<ul style="list-style-type: none"> • 350lbs 	<ul style="list-style-type: none"> • W 	TLFM680
	<p>Titan Swing Bracket</p> <p>The TLFM685 pairs with the 72" or 36" Horizontal Rail (TLSG620, TLSG625), allowing the Link System to terminate at a wall.</p>	<ul style="list-style-type: none"> • 350lbs 	<ul style="list-style-type: none"> • W 	TLFM685

Technical Specifications & Accessories