




# SG6101 SafetyGrip® GRAB POLE WITH ROTOBAR

Page 1 of 2



**THE HEIGHT OF THE ROTOBAR TOP BAR FROM THE FLOOR IS FULLY ADJUSTABLE FROM 27" to 43"**

**THE CEILING HEIGHT OF THIS GRAB POLE IS ADJUSTABLE FROM 7'-9" to 8'-3"**

**THE CEILING HEIGHT OF THIS GRAB POLE CAN BE INCREASED UP TO 10'-3" WITH OUR SG604 HEIGHT EXTENSION KIT**


**WEIGHT RATING: 350lbs**

**NOTE: This product is not intended to take the full weight of the user**



# SG6101 SafetyGrip® GRAB POLE WITH ROTOBAR

Page 2 of 2



**SG6101W  
WHITE**



**SG6101G  
GREY**



**SG6101B  
BLACK**